

# BD Introsyte™ Autoguard™ Shielded Introducer

*Preferred Performance*  
*Proven Needlestick Protection*



Indispensable to  
human health



# BD Introsyte™ Autoguard™ Shielded Introducer

## For PICC and Midline Placement

The BD Introsyte Autoguard shielded introducer incorporates the BD Autoguard needle shielding technology, proven to reduce IV stylet-related needlestick injuries by 95 percent.\* A simple push of the button safely retracts the needle, with little change in user technique. The BD Introsyte Autoguard shielded introducer provides you with a safe and effective introduction method for successful PICC and midline placement. In addition, the BD Introsyte-N™ Autoguard™ shielded introducer includes a notched needle for earlier visualization of blood return and is specifically designed for the neonatal and pediatric patient population.



## Superior Insertion Features

The BD Introsyte Autoguard family of introducers incorporates superior insertion characteristics for improved access. BD's high quality cannula tipping process creates a minimal lie distance between the needle bevel and the sheath. This proprietary tipping technology with superior needle lubrication enhances first-stick proficiency and ease of insertion.



BD Autoguard™ technology is the **only** needle-shielding technology that has been clinically proven to make a *"...marked and significant reduction in IV stylet-related injuries..."*\*

# BD Extended Dwell Catheters – Quality and Reliability You Can Count On

## BD First PICC™ and BD First Midcath™ Extended Dwell Catheters

BD First PICC and BD First Midcath are high-quality silicone catheters designed for easy and safe placement.

- Available in a wide range of single- and dual-lumen gauge sizes
- Soft silicone integrated extension set enhances patient comfort and provides increased mobility
- Protective sleeve prevents accidental contamination of the catheter
- Pre-attached T-port extension allows you to flush the catheter before and after insertion and provides added security against stylet movement
- Numbered catheter markings allow for easy measurement and catheter trimming prior to insertion
- Latex- and DEHP-free

### Catheter Components:

- 1 BD First PICC or BD First Midcath silicone catheter with hydrophilic-coated stylet
- 1 T-port extension
- 1 Tape measure
- 1 Patient information booklet
- 1 Patient information card
- 1 Instructions for use



## BD L-Cath™ Extended Dwell Catheters

BD L-Cath PICC and midline catheters are designed for easy and safe placement. The polyurethane catheter material provides strength, versatility and increased flow rates.

- Available in a wide range of single- and dual-lumen gauge sizes
- Protective sleeve prevents accidental contamination of the catheter
- Pre-attached Y-site allows you to flush the catheter before and after insertion and provides added security against stylet movement
- Numbered catheter markings allow for easy measurement and catheter trimming prior to insertion
- Flexible, braided, stainless steel wire stylet facilitates easy catheter insertion
- Trim tool ensures clean and accurate trimming of the catheter to desired length
- Latex- and DEHP-free

### Catheter Components:

- 1 BD L-Cath PICC or midline polyurethane catheter with braided stainless steel wire stylet
- 1 Trim tool\*\*
- 1 Instructions for use



\* Mendelson M.H., Chen L.B.Y., Finkelstein L.E., Bailey E., Kogan G., Mt. Sinai Medical Center, New York, NY. "Preliminary Evaluation of a Safety IV Catheter (Insyte Autoguard, Becton, Dickinson) Using the Centers for Disease Control and Prevention (CDC) National Surveillance System for Hospital Healthcare Workers Database." Abstr. Centers for Disease Control and Prevention - 4th Decennial International Conference on Nosocomial and Healthcare-Associated Infections. March 2000.

\*\* Not included with midline catheters

BD Introsyte™ Autoguard™ Shielded Introducers		
384090	24 Ga/1.9F Introsyte-N™ Autoguard™ Introducer	10 ea/cs
384091	20 Ga/3F Introsyte Autoguard Introducer	10 ea/cs
384092	18 Ga/4F Introsyte Autoguard Introducer	10 ea/cs
384093	16 Ga/5F Introsyte Autoguard Introducer	10 ea/cs

BD First PICC™ Silicone Basic Kits		
384115	20 Ga x 65 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs
384171	18 Ga x 65 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs
384177	16 Ga x 65 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs
384183	16 Ga x 65 cm D/L Basic w/ Introsyte Autoguard	5 ea/cs

BD L-Cath™ PICC Polyurethane Basic Kits		
384566	20 Ga x 60 cm D/L Basic w/ Introsyte Autoguard	5 ea/cs
384481	3F x 60 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs
384490	18 Ga x 60 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs
384575	18 Ga x 60 cm D/L Basic w/ Introsyte Autoguard	5 ea/cs
384499	16 Ga x 60 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs
384584	16 Ga x 60 cm D/L Basic w/ Introsyte Autoguard	5 ea/cs

BD First Midcath™ Silicone Basic Kits		
384304	20 Ga x 20 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs
384308	18 Ga x 20 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs

BD L-Cath™ Midline Polyurethane Basic Kits		
384631	3F x 20 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs
384635	18 Ga x 20 cm S/L Basic w/ Introsyte Autoguard	5 ea/cs

BD First PICC™ Silicone Intermediate Procedure Kits		
384117	20 Ga x 65 cm S/L Intermediate w/ Introsyte Autoguard	5 ea/cs
384173	18 Ga x 65 cm S/L Intermediate w/ Introsyte Autoguard	5 ea/cs
384178	16 Ga x 65 cm S/L Intermediate w/ Introsyte Autoguard	5 ea/cs
384185	16 Ga x 65 cm D/L Intermediate w/ Introsyte Autoguard	5 ea/cs

BD L-Cath™ PICC Polyurethane Complete Procedure Kits		
384572	20 Ga x 60 cm D/L Complete w/ Introsyte Autoguard	5 ea/cs
384487	3F x 60 cm S/L Complete w/ Introsyte Autoguard	5 ea/cs
384496	18 Ga x 60 cm S/L Complete w/ Introsyte Autoguard	5 ea/cs
384581	18 Ga x 60 cm D/L Complete w/ Introsyte Autoguard	5 ea/cs
384563	16 Ga x 60 cm S/L Complete w/ Introsyte Autoguard	5 ea/cs
384590	16 Ga x 60 cm D/L Complete w/ Introsyte Autoguard	5 ea/cs

BD First Midcath™ Silicone Intermediate Procedure Kits		
384306	20 Ga x 20 cm S/L Intermediate w/ Introsyte Autoguard	5 ea/cs
384310	18 Ga x 20 cm S/L Intermediate w/ Introsyte Autoguard	5 ea/cs

BD L-Cath™ Midline Polyurethane Intermediate Procedure Kits		
384633	3F x 20 cm S/L Intermediate w/ Introsyte Autoguard	5 ea/cs
384637	18 Ga x 20 cm S/L Intermediate w/ Introsyte Autoguard	5 ea/cs

**For More Information on  
BD Extended Dwell Catheter Products:  
Call 888-237-2762**



**BD Medical Systems**  
9450 South State Street  
Sandy, UT 84070  
www.bd.com

**Catheter Only**

First PICC, Silicone Catheter Only		
384232	24 Ga/1.9F S/L Silicone PICC, 50 cm	5 ea/cs
384156	20 Ga/3F S/L Silicone PICC, 65 cm	5 ea/cs
384157	18 Ga/4F S/L Silicone PICC, 65 cm	5 ea/cs
384158	16 Ga/5F S/L Silicone PICC, 65 cm	5 ea/cs
384159	16 Ga/5F D/L Silicone PICC, 65 cm	5 ea/cs

L-Cath PICC, Polyurethane Catheter Only		
384539	24 Ga/1.9F S/L Poly PICC, 30 cm	10 ea/cs
384463	20 Ga/3F S/L Poly PICC, 60 cm	10 ea/cs
384464	18 Ga/3.5F S/L Poly PICC, 60 cm	10 ea/cs
384465	16 Ga/5F S/L Poly PICC, 60 cm	10 ea/cs
384466	20 Ga/2.6F D/L Poly PICC, 60 cm	10 ea/cs
384467	18 Ga/3.5F D/L Poly PICC, 60 cm	10 ea/cs
384468	16 Ga/5F D/L Poly PICC, 60 cm	10 ea/cs

First Midcath, Silicone Catheter Only		
384357	20 Ga/3F S/L Silicone Midline, 20 cm	5 ea/cs
384358	18 Ga/4F S/L Silicone Midline, 20 cm	5 ea/cs

L-Cath Midline, Polyurethane Catheter Only		
384620	20 Ga/3F S/L Poly Midline, 20 cm	10 ea/cs
384621	18 Ga/3.5F S/L Poly Midline, 20 cm	10 ea/cs

**BD Procedure Trays**

384825	Neonatal Procedure Tray for PICC/midline insertion	5 ea/cs
384826	Intermediate Procedure Tray for PICC/midline insertion	5 ea/cs
384827	Complete Procedure Tray for PICC/midline insertion	5 ea/cs

384825 BD Neonatal Procedure Kit Contents:	384826 BD Intermediate Procedure Tray Contents:	384827 BD Complete Procedure Tray Contents:
1 Face mask with eye shield	2 Isolation masks	1 Face mask with eye shield
1 Tape measure	1 Tape measure	1 Isolation mask
1 Poly-lined drapes	1 Poly-lined drape	1 Tape measure
1 Fenestrated drape	1 Fenestrated drape	1 Poly-lined drape
1 Iodophor PVP swabsticks (3 per pack)	1 Iodophor PVP swabsticks (3 per pack)	1 Fenestrated drape
1 Alcohol swabsticks (3 per pack)	1 Alcohol swabsticks (3 per pack)	1 Iodophor PVP swabsticks (3 per pack)
1 Sterile normal saline wipe	1 Sterile normal saline wipe	1 Alcohol swabsticks (3 per pack)
1 Tourniquet	1 Tourniquet	1 Sterile normal saline wipe
2 4x4 Absorbent pads	4 4x4 Absorbent pads	1 Tourniquet
2 2x2 Absorbent pads	2 2x2 Absorbent pads	2 Gloves
1 BD™ Syringe, 10 cc	2 BD™ Syringes, 10 cc	1 BD SafetyGlide™ syringe with needle, 1 cc
1 BD Posiflow™ Multi-dose vial adapter	1 BD Posiflow™ Multi-dose vial adapter	4 4x4 Absorbent pads
1 BD Posiflow™ luer-activated valve	1 BD Posiflow™ luer-activated valve	2 2x2 Absorbent pads
1 Scissors	1 Scissors	2 BD™ Syringes, 10 cc
1 Forceps, with tip protection	1 Forceps, with tip protection	2 Ampules, 10 ml sodium chloride 0.9%
1 Extension set	4 Foam tape strips	1 Filter straw
3 Foam tape strips	1 Roll Transpore™ tape	1 BD Posiflow™ luer-activated valve
1 Roll Transpore™ tape	1 Skin protectant swabstick	1 Scissors
1 Skin protectant swabstick	1 Tegaderm™ HP transparent dressing	1 Forceps, with tip protection
1 Tegaderm™ HP transparent dressing		1 Extension set
		1 StatLock® catheter securement
		1 First StatLock® catheter securement
		4 Foam tape strips
		1 Roll Transpore™ tape
		1 Skin protectant swabstick
		1 Tegaderm™ HP transparent dressing

StatLock is a trademark of Venetec International, Inc.  
Tegaderm and Transpore are trademarks of 3M.  
BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. ©2003 BD