

# MONOJECT® Safety Syringe



## Safety By Design

Complies with OSHA  
Bloodborne Pathogens Standard  
for Engineering Controls

## Now Available in SoftPack

Tyco Healthcare manufactures only the highest quality products that provide for your safety. The MONOJECT Safety Syringe is proven effective and preferred by its millions of users each day. Our exclusive design helps keep clinicians safe with maximum clinical flexibility.

### Secure Safety Shield

- The robust safety shield is strong and secure, covering the needle completely during transport and after usage.

### Transport Position

- MONOJECT's exclusive transport position allows clinicians to safely transport medication and eliminates the need to recap the needle. Only the MONOJECT Safety Syringe has this feature.

### Intermittent Drug Delivery

- The exclusive transport position lets clinicians temporarily cover the needle with the shield instead of recapping it while titrating medication, keeping fingers away from the needle.

### Shield Activation with Fluid in the Barrel

- The MONOJECT safety shield can be activated with fluid still in the barrel and with no risk of splatter.


### Full Product Offering

- Safety needs do not come in just one size. That is why MONOJECT offers the 3/10cc and 1/2cc Safety Insulin in addition to 1cc, 3cc, 6cc, and 12cc sizes.

**tyco**  
Healthcare

**Kendall**

**SharpSafety**  
Needles & Syringes • Sharps Disposal • Blood Collection



# MONOJECT® Safety Syringe



## MONOJECT® Safety Syringe

Item Code	Description	Needle Gauge & Length	Box/Case Quantity
8881511201	1cc TB	28 GA. x 1/2"	100/500
8881511235	1cc TB	25 GA. x 5/8"	100/500
8881511110	1cc Insulin	29 GA. x 1/2"	100/500
8881511136	1/2cc Insulin	29 GA. x 1/2"	100/500
8881511144	3/10cc Insulin	29 GA. x 1/2"	100/500
8881533908	3cc Safety Syringe Only		100/500
8881533056	3cc	20 GA. x 1 1/2"	100/500
8881533130	3cc	21 GA. x 1"	100/500
8881533155	3cc	21 GA. x 1 1/2"	100/500
8881533239	3cc	22 GA. x 1"	100/500
8881533254	3cc	22 GA. x 1 1/2"	100/500
8881533338	3cc	23 GA. x 1"	100/500
8881533510	3cc	25 GA. x 5/8"	100/500
8881566007	6cc Safety Syringe Only		50/300
8881566015	6cc	20 GA. x 1 1/2"	50/300
8881566023	6cc	21 GA. x 1 1/2"	50/300
8881522000	12cc Safety Syringe Only		50/300
8881522018	12cc	20 GA. x 1 1/2"	50/300
8881522026	12cc	21GA. x 1 1/2"	50/300
<i>Now Available In SoftPack</i>			
8881033908	3cc Safety Syringe Only		80/480
8881033056	3cc	20 GA. x 1 1/2"	80/480
8881033130	3cc	21 GA. x 1"	80/480
8881033155	3cc	21 GA. x 1 1/2"	80/480
8881033239	3cc	22 GA. x 1"	80/480
8881033254	3cc	22 GA. x 1 1/2"	80/480
8881033338	3cc	23 GA. x 1"	80/480
8881533510	3cc	25 GA. x 5/8"	80/480

MONOJECT® is a trademark of Sherwood Services AG.

**tyco**  
Healthcare

**Kendall**

15 Hampshire Street  
Mansfield, MA 02048  
1-800-962-9888  
508-261-8000

[www.kendallhq.com](http://www.kendallhq.com)

For International locations please see  
[www.tycohealthcare.com](http://www.tycohealthcare.com)

H-5315 5M 10/01