





Proxima

Surgical Drapes

Designed,
Tested and
Trusted
by O.R
Professionals



Proxima Surgical Drapes



Better by design

When we developed our Proxima line of surgical drapes we went straight to the people who use them. We spoke to O.R. professionals across the country to get a sense of what features they liked and disliked about other surgical drapes and then worked to develop Proxima — an innovative line of surgical drapes that you can use with confidence.

The drapeability of Softesse™ medical fabric from DuPont™

Proxima drapes feature Softesse medical fabric from DuPont, an industry standard for more than 25 years. The fabric gives Proxima drapes their soft, cloth-like feel that conforms easily to the body.

Velcro-style line holders

Proxima drapes feature velcro-style line holders to place surgical lines and help keep them clear of the fenestration. The tabs are offset to make lifting them easy, even when wearing gloves. This is one of the key features suggested to us by clinicians across the country when we developed the Proxima line of drapes.

Premium tape and incise film

We use premium tape and incise film exclusively on all Proxima drapes. Our tape is manually placed flush to the fenestration to help prevent pooling of blood and fluid underneath the drape. Because the tape is placed manually, it allows the tape to be closer to the fenestration compared to other drapes.

Larger impervious reinforcement zone

Proxima fenestrations and splits are surrounded with a large impervious reinforcement. This area absorbs fluids and is plastic-backed to make the area impervious to blood and body fluids as well as bacterial transmission. All Proxima reinforced zones pass ASTM F1670 and F1671.

Wider anesthesia screen

The anesthesia screen on Proxima drapes is 102" wide, allowing you to drape a 6'9" person from fingertip to fingertip. This virtually eliminates the need to pull half-sheets or three-quarter sheets to create a sterile field. Over time, this feature can be a real money-saver.

Instrument pad for traction control

The instrument pad on Proxima drapes is a section of open-celled foam that provides a convenient, non-slip area for placing instruments during surgery. The pad is also plastic-backed, making it impervious to fluid transfer.

Medline's Proxima surgical drape line is the broadest and most extensive in the industry.

The Proxima “Mini-Drape” is a representation of a Proxima laparotomy drape. It shows several key features including the velcro-style line holders, premium tape, our larger impervious reinforcement zone, the instrument pad, and the Softesse™ medical fabric from DuPont.™



Marianne Finlay
Director of Surgical Services
Vista Health

“When Medline said they could save us money, I listened. But when I experienced the quality of Proxima gowns and drapes, I believed.”



“Whenever we introduce a product to the O.R., there’s always a concern that perhaps we’re sacrificing quality, especially if we’re saving money.

I’m pleased to report that when we started using Proxima drapes and gowns from Medline, we not only saved a great deal of money, we upgraded the quality level from our previous line of drapes and gowns.

It’s a win-win situation.”

Ask your Medline representative to see our Proxima catalog — featuring the most extensive line of drapes and gowns in the industry.